→ BPL_TEST2_Fedbatch script with PyFMI ver 2.9.8

The key library PyFMI ver 2.9.8 is installed.

After the installation a small application BPL_CHO_Fedbatch is loaded and run. You can continue with this example if you like.

```
!lsb release -a # Actual VM Ubuntu version used by Google
    No LSB modules are available.
    Distributor ID: Ubuntu
    Description:
                    Ubuntu 20.04.5 LTS
    Release:
                    20.04
    Codename:
                    focal
!python --version
    Pvt.hon 3.9.16
%env PYTHONPATH=
    env: PYTHONPATH=
!wget https://repo.anaconda.com/miniconda/Miniconda3-py39_23.1.0-1-Linux-x86_64.sh
!chmod +x Miniconda3-py39_23.1.0-1-Linux-x86_64.sh
!bash ./Miniconda3-py39_23.1.0-1-Linux-x86_64.sh -b -f -p /usr/local
import sys
sys.path.append('/usr/local/lib/python3.9/site-packages/')
    --2023-03-29 11:34:03-- https://repo.anaconda.com/miniconda/Miniconda3-py39_23.1.0-1-Linux-x86_64.sh
    Resolving repo.anaconda.com (repo.anaconda.com)... 104.16.130.3, 104.16.131.3, 2606:4700::6810:8203, ...
    Connecting to repo.anaconda.com (repo.anaconda.com) | 104.16.130.3 | :443... connected.
    HTTP request sent, awaiting response... 200 OK
    Length: 69888122 (67M) [application/x-sh]
    Saving to: 'Miniconda3-py39 23.1.0-1-Linux-x86 64.sh'
    Miniconda3-py39 23. 100%[===========] 66.65M
                                                             187MB/s
                                                                         in 0.4s
    2023-03-29 11:34:04 (187 MB/s) - 'Miniconda3-py39_23.1.0-1-Linux-x86_64.sh' saved [69888122/69888122]
    PREFIX=/usr/local
    Unpacking payload ...
    Installing base environment...
    Downloading and Extracting Packages
    Downloading and Extracting Packages
    Preparing transaction: done
    Executing transaction: done
    installation finished.
!conda update -n base -c defaults conda --yes
```

```
Preparing transaction: done
Verifying transaction: done
Executing transaction: done
!conda --version
!python --version
```

conda 23.3.0 Python 3.9.16

!conda install -c conda-forge pyfmi --yes # Install the key package

```
Preparing transaction: done
Verifying transaction: done
Executing transaction: done
```


→ BPL_TEST2_Fedbatch setup

Now specific installation and the run simulations. Start with connecting to Github. Then upload the two files:

- FMU BPL_TEST2_Fedbatch_linux_om_me.fmu
- Setup-file BPL_TEST2_Fedbatch_explore.me.py

```
%%bash
git clone https://github.com/janpeter19/BPL_TEST2_Fedbatch
    Cloning into 'BPL TEST2 Fedbatch'...
%cd BPL_TEST2_Fedbatch
    /content/BPL_TEST2_Fedbatch
run -i BPL_TEST2_Fedbatch_explore.py
    Linux - run FMU pre-comiled OpenModelica 1.21.0
    Model for bioreactor has been setup. Key commands:
     - par()
                  - change of parameters and initial values
     - init()
                   - change initial values only
                   - simulate and plot
     - simu()
                  - make a new plot
     - newplot()
                   - show plot from previous simulation
     - show()
                   - display parameters and initial values from the last simulation
     - disp()
     - describe() - describe culture, broth, parameters, variables with values / units
    Note that both disp() and describe() takes values from the last simulation
    Brief information about a command by help(), eg help(simu)
    Key system information is listed with the command system_info()
%matplotlib inline
plt.rcParams['figure.figsize'] = [25/2.54, 20/2.54]
import warnings
warnings.filterwarnings("ignore")
```


▼ BPL TEST2 Fedbatch - demo



disp(mode='long')

```
bioreactor.V_0: V_0: 1.0
bioreactor.m_0[1]: VX_0: 1.0
bioreactor.m_0[2]: VS_0: 10.0
bioreactor.culture.Y: Y: 0.5
bioreactor.culture.GSmax: gSmax: 1.0
bioreactor.culture.Ks: Ks: 0.1
feedtank.C_in[2]: feedtank.S_in: 300.0
feedtank.V_0: feedtank.V_0: 10.0
dosagescheme.mu_feed: mu_feed: 0.1
dosagescheme.t_start: t_start: 3.0
dosagescheme.F_start: F_start: 0.001
dosagescheme.F_max: F_max: 0.3
```

A more typical feed scheme for the culture at hand newplot(plotType='TimeSeries') par(t_start=4, F_start=0.008, mu_feed=0.2, F_max=0.1) simu(20)



```
disp('culture')

Draw Y: 0.5
    qSmax: 1.0
    Ks: 0.1

describe('mu')

    Cell specific growth rate variable: 0.12 [ 1/h ]

describe('parts')
    ['bioreactor', 'bioreactor.culture', 'dosagescheme', 'feedtank']

describe('MSL')
    MSL: 3.2.3 - used components: RealInput, RealOutput

system_info()
```

system_IIIIO

```
System information
-OS: Linux
-Python: 3.9.16
-Scipy: not installed in the notebook
-PyFMI: 2.10.0
-FMU by: OpenModelica Compiler OpenModelica 1.21.0~dev-185-g9d983b8
-FMI: 2.0
-Type: FMUModelME2
-Name: BPL_TEST2.Fedbatch
-Generated: 2023-01-19T09:29:14Z
-MSL: 3.2.3
-Description: Bioprocess Library version 2.1.1-beta
-Interaction: FMU-explore version 0.9.7
```

Colab paid products - Cancel contracts here

✓ 0s completed at 13:37