## BPL\_TEST2\_Fedbatch script with FMPy

The key library FMPy is installed.

openssl

After the installation a small application BPL\_TEST2\_Fedbatch is loaded and run. You can continue with this example if you like.

```
!lsb_release -a # Actual VM Ubuntu version used by Google

→ No LSB modules are available.
    Distributor ID: Ubuntu
    Description:
                     Ubuntu 22.04.3 LTS
    Release:
                     22.04
    Codename:
                     jammy
%env PYTHONPATH=
→ env: PYTHONPATH=
!wget https://repo.anaconda.com/miniconda/Miniconda3-py312_24.3.0-0-Linux-x86_64.sh
!chmod +x Miniconda3-py312_24.3.0-0-Linux-x86_64.sh
!bash ./Miniconda3-py312_24.3.0-0-Linux-x86_64.sh -b -f -p /usr/local
import sys
sys.path.append('/usr/local/lib/python3.12/site-packages/')
   --2024-08-13 08:15:18-- https://repo.anaconda.com/miniconda/Miniconda3-py312 24.3.0-0-Linux-x86 64.sh
    Resolving repo.anaconda.com (repo.anaconda.com)... 104.16.191.158, 104.16.32.241, 2606:4700::6810:bf9e, ...
    Connecting to repo.anaconda.com (repo.anaconda.com)|104.16.191.158|:443... connected.
    HTTP request sent, awaiting response... 200 OK
Length: 143351488 (137M) [application/octet-stream]
    Saving to: 'Miniconda3-py312_24.3.0-0-Linux-x86_64.sh'
    Miniconda3-py312_24 100%[==========] 136.71M 83.3MB/s
    2024-08-13 08:15:20 (83.3 MB/s) - 'Miniconda3-py312_24.3.0-0-Linux-x86_64.sh' saved [143351488/143351488]
    PREFIX=/usr/local
    Unpacking payload ...
    Installing base environment...
    Preparing transaction: ...working... done
    Executing transaction: ...working... done
    installation finished.
!conda update -n base -c defaults conda --yes
\rightarrow
                                              2024.3.11-h06a4308_0 --> 2024.7.2-h06a4308_0
      ca-certificates
                                          2024.2.2-py312h06a4308_0 --> 2024.7.4-py312h06a4308_0
       certifi
                                            24.3.0-py312h06a4308_0 --> 24.7.1-py312h06a4308_0
3.0.13-h7f8727e_0 --> 3.0.14-h5eee18b_0
       conda
```

CUIIUd-24./.1 | 1.2 PID | : 100% 1.0/1 [000:000, 1./01./5]

frozendict-2.4.2 | 36 KB | : 100% 1.0/1 [00:00<00:00, 1.69it/s]

Preparing transaction: done Verifying transaction: done Executing transaction: done

!conda --version
!python --version

conda 24.7.1 Python 3.12.2

!conda install -c conda-forge fmpy --yes # Install the key package

 $\overline{\mathbf{T}}$ 

Preparing transaction: done Verifying transaction: done Executing transaction: done

```
#!conda install matplotlib --yes
#!conda install scipy --yes
#!conda install xlrd --yes
#!conda install openpyxl --yes
```

## BPL\_TEST2\_Fedbatch setup

Now specific installation and the run simulations. Start with connecting to Github. Then upload the two files:

```
• FMU - BPL_TEST2_Fedbatch_linux_om_me.fmu
```

· Setup-file - BPL\_TEST2\_Fedbatch\_fmpy\_explore.py

```
git clone https://github.com/janpeter19/BPL_TEST2_Fedbatch
→ Cloning into 'BPL_TEST2_Fedbatch'...
%cd BPL_TEST2_Fedbatch
→ /content/BPL_TEST2_Fedbatch
run -i BPL_TEST2_Fedbatch_fmpy_explore.py
→ Linux - run FMU pre-comiled OpenModelica 1.23.0-dev
    Model for bioreactor has been setup. Key commands:
                   - change of parameters and initial values
     - par()
     - init()
                   - change initial values only
     - simu()

    simulate and plot

    make a new plot


     - newplot()

    show plot from previous simulation


      - show()
                   - display parameters and initial values from the last simulation
     - disp()
     - describe() - describe culture, broth, parameters, variables with values/units
    Note that both disp() and describe() takes values from the last simulation
    and the command process_diagram() brings up the main configuration
    Brief information about a command by help(), eg help(simu)
    Key system information is listed with the command system_info()
%matplotlib inline
plt.rcParams['figure.figsize'] = [25/2.54, 20/2.54]
import warnings
warnings.filterwarnings("ignore")
BPL_TEST2_Fedbatch - demo
describe('culture'); print(); #describe('liquidphase')
                                                                                        # Pump schedule parameter

→ Simplified text book model – only substrate S and cell concentration X
process_diagram()
```

→ No processDiagram.png file in the FMU, but try the file on disk.




# Simulation with default values of the process
newplot(plotType='TimeSeries')
simu(20)



disp(mode='long')

```
bioreactor.V_start : V_start : 1.0
bioreactor.m_start[1] : VX_start : 0.0
bioreactor.m_start[2] : VS_start : 0.0
bioreactor.culture.Y : Y : 0.5
bioreactor.culture.qSmax : qSmax : 1.0
bioreactor.culture.Ks : Ks : 0.1
feedtank.c_in[2] : feedtank.S_in : 0.0
feedtank.V_start : feedtank.V_start : 100.0
dosagescheme.mu_feed : mu_feed : 0.2
dosagescheme.t_startExp : t_startExp : 2.0
dosagescheme.F_startExp : F_startExp : 0.12
dosagescheme.F_max : F_max : 3.0
# A more typical feed scheme for the culture at hand
newplot(plotType='TimeSeries')
par(t_startExp=4, F_startExp=0.008, mu_feed=0.2, F_max=0.1)
simu(20)
```





```
disp('culture')
   Y: 0.5
    qSmax : 1.0
    Ks : 0.1
describe('mu')
→ Cell specific growth rate variable : 0.12 [ 1/h ]
describe('parts')
describe('MSL')

→ MSL: 3.2.3 – used components: RealInput, RealOutput
system_info()
    System information
     -OS: Linux
     -Python: 3.10.12
     -Scipy: not installed in the notebook
     -FMPy: 0.3.20
     -FMU by: OpenModelica Compiler OpenModelica 1.23.0~dev.beta.1-1-g379f714
     -FMI: 2.0
     -Type: ME
-Name: BPL_TEST2.Fedbatch
     -Generated: 2024-05-20T13:38:03Z
     -MSL: 3.2.3
     -Description: Bioprocess Library version 2.2.0
     -Interaction: FMU-explore for FMPy version 1.0.1
Start coding or generate with AI.
```